

Frecuencia de Mutación

Manual de Usuario

Sección II.5: Frecuencia de Mutación (p105 - p109)

Consulta el manual completo [aquí](#).

Archivo de Ejemplo

Nombre del archivo: FRECUENCIA DE MUTACIÓN RT VIH-1.csv

Selección de Parámetros

- **Nucleótidos/Aminoácidos:**
 - *Aminoácidos*
- **Alinear:**
 - *No marcar*
- **Referencia:**
 - *PISPIETVPVKLPGMDGPKVKQWPLTEEKIKALVEICTEMEKEGKISKIGPENPYNTPVFAIKKK
DSTKWRKLVDFRELNKRTQDFWEVQLGIPHPAGLKKKSVTVLDVGDAYFSVPLDEDFRKYTAF
TIPSINNETPGIRYQYNVLPQGWKGSPAIFQSSMTKILEPFRKQNPDIVIYQYMDDLIVGSDLEIG
QHRTKIEELRQHLLRWGLTTPDKKHQKEPPFLWMGYELHPDKWTVQPIVLPEKDSWTVNDIQ
KLVGKLNWASQIYPGIKVRQLCKLLRGTKALTEVIPLTEEALELAENREILKEPVHGVYDPSKDLI
AEIQKQGQGQWTYQIYQEPFKNLKTGKYARMRGAHTNDVKQLTEAVQKITTESIVIWGKTPKF
KLPIQKETWETWWTEYWQATWIPEWEFVNTPLVLKWLWYQLEKEPIVGAETF*



The screenshot shows the EpiMolBio v.0.1 web application interface. The title bar indicates the version. The main header is 'FRECUENCIA DE MUTACIÓN' with a DNA helix icon. A left sidebar contains navigation links: EpiMolBio, INICIO, VIRUS, VARIABILIDAD (highlighted), HOMOLOGÍA, RASTREADOR, ALINEAMIENTOS, and HERRAMIENTAS. At the bottom of the sidebar are flags for Spain and the UK. The main content area has a green header. Below it, there's an 'Entrada' (Input) section with a text box and a '< Atrás' button. The 'Referencia' (Reference) section has a text box containing 'WEFVNTPLVLKWLWYQLEKEPIVGAETF'. Below the input and reference sections are checkboxes for 'Nucleótidos' (unchecked), 'Aminoácidos' (checked), and 'Alinear' (unchecked). There's a 'Salida' (Output) section with a text box. At the bottom, there are icons for information and help, and two buttons: 'Calcular' (Calculate) and 'Cancelar' (Cancel).