

## Alignments Delete Insertions

**It allows automatically removing insertions from a sequence with respect to a reference with gaps after performing the alignment.**

- ① Select the input folder. It should only contain .fasta files of aligned NT or AA sequences to be analyzed.
- ② Enter the reference sequence with gaps and without line breaks or spaces.
- ③ Select the output folder without entering the file name. It will be automatically named as 'Insertions\_Deleted\_InputFileName.fasta'.
- ④ Click Calculate.

A .fasta file will be generated with the sequences having insertions removed.

More information about this function can be found in the User Manual, section:

V. ALIGNMENTS V.3 DELETE INSERTIONS