

Manual de Usuario

Sección I.1.A: Mutaciones de Resistencia - Transmitidas (p9 - p11, p14 - p18)


Consulta el manual completo [aquí](#).

Archivo de Ejemplo

Nombre del archivo: SDRM ITIAN VIH-1 LISTA.html

Selección de Parámetros

- **Tipo de SDRM:**
 - *SDRM - ITIAN*
- **Formato de Salida:**
 - *Lista*



The screenshot shows the EpiMolBio v.0.1 web application interface. The left sidebar contains a menu with the following items: EpiMolBio, INICIO, VIRUS (highlighted in green), VARIABILIDAD, HOMOLOGÍA, RASTREADOR, ALINEAMIENTOS, and HERRAMIENTAS. At the bottom of the sidebar are flags for Spain and the United Kingdom. The main content area is titled 'TRANSMITIDAS' and features a green header bar with a logo. Below the header, there are two input fields: 'Entrada' and 'Salida', both with the text 'Entrada' and 'Salida' respectively. To the right of the 'Entrada' field is a button labeled '< Atrás'. Below these fields are two dropdown menus: 'Tipo de SDRM' (set to 'SDRM-ITIAN') and 'Formato de Salida' (set to 'Lista'). At the bottom of the main area are two buttons: 'Calcular' and 'Cancelar'. There are also two small icons (an information icon and a warning icon) located below the 'Salida' field.